

Features and functions of the Cardinal Health™ Kangaroo OMNI™ Enteral Feeding Pump

The Cardinal Health™ Kangaroo OMNI™ Enteral Feeding Pump is designed to meet feeding needs from hospital to home. The Kangaroo OMNI™ pump is small, light and easy to transport. It offers the automated feed/flush settings the Kangaroo™ Portfolio is known for, with the added benefit of delivering thick, homogenous, blenderized formulas categorized as Level 2, 3 or 4 drinks within the IDDSI framework. IDDSI provides standardized definitions to describe fluid thickness.†

See how the pump's flexibility and user-friendly features support your enteral feeding needs from hospital to home.

1 Automated feed/flush technology: Automated water flushes deliver hydration and help avoid clogs in the feeding tube.

2 Timed pause: You can hold feeding for a specified amount of time from 5 to 240 minutes without turning off the pump or re-entering pump settings.

3 Lock screen: The lock screen disables keypad buttons to help prevent enabling pump features during a feeding.

4 Feeding history*: The pump displays and stores up to 30 days of feeding history, the longest of any Kangaroo™ Enteral Feeding Pump.

5 Accuracy: The Kangaroo OMNI™ Enteral Feeding Pump has an accuracy of +/- 5% or 0.5 mL when used with standard feed only or feed/flush sets and +/- 10% or 1.0 mL when used with thick formula feeding sets or flushing.†

6 Night mode: Users can dim the screen brightness in night mode to help reduce disruptions while sleeping. The Kangaroo OMNI™ Pump is the quietest Kangaroo™ Feed and Flush Pump.

7 LED indicator lights: Red, yellow and green LED lights visually indicate the pump status.

10 Lithium-ion rechargeable battery: The longest battery life of any Kangaroo™ Product, with a battery life up to 20 hours.*

11 Portability: Designed to move with you, the Kangaroo OMNI™ Pump has a rubberized handle, weighs 1.2 lbs (550 g) and can operate in any orientation, even in the Kangaroo™ Backpack.

12 Cassette design: The Kangaroo OMNI™ Feeding Sets use a color-coding system to easily identify each cassette. The user can load the Kangaroo OMNI™ Feeding Set with one hand.

13 Intuitive display: A large, color LCD display features step-by-step prompts, animated illustrations to guide pump operation and key notifications to improve user experience.

8 Interruption monitor screen: Tracks the delivered amount of feed and hydration and will calculate and display the amount of missed nutrition if the pump is turned off. This allows for an objective approach to assessing prescription adherence. This notice screen will be displayed if the pump has been turned off, stopped or flushing for longer than 30 minutes.

9 Thick formula delivery: The Kangaroo OMNI™ can deliver thick or blenderized formula — categorized as Level 2, 3 or 4 drinks within the IDDSI framework — which provides a framework for drink textures.



To learn more contact your Cardinal Health sales representative or visit cardinalhealth.com/kangaroo

Important information: Prior to use, refer to the Instructions For Use supplied with this device for indications for use, suggested procedure, warnings and cautions.

*when running at 125 mL/hr

†Thick formula is defined as enteral fluids of smooth consistency that would be categorized as level 2, 3 or 4 drinks within the IDDSI framework.

‡This feature is not being displayed in this image.

© 2024 Cardinal Health. All Rights Reserved. CARDINAL HEALTH, the Cardinal Health LOGO, KANGAROO, the Kangaroo LOGO and Kangaroo OMNI are trademarks of Cardinal Health and may be registered in the US and/or in other countries. All other marks are the property of their respective owners. Lit. No. 2MS23-2518268-06 (01/2024)

